

HTS2

HTS2 - Online VT Pursuits Doctor Manual

© Copyright HTS Inc. 2019



HTS2 - Online VT

Developed By

Dr. Jeffery Cooper, OD, FCOVD
Rodney K. Bortel
Ted L. Conner

Pursuits

When in the default Daily Therapy Protocol, Pursuits will run for 3 minutes each session.

TODAY'S ASSIGNMENT	MIN.
Pursuits	03:00
Saccades	03:00
Divergence	07:00
Convergence	07:00
Accommodative Rock	05:00

If all of the Stars have been earned for Pursuits and you wish your patient to continue Pursuits, you may place the Exercise in Maintenance Mode. A Exercise in Maintenance Mode will always be displayed in the patients Daily Therapy Protocol until removed from Maintenance. Maintenance Mode may be accessed by Clicking Existing Patients from the Navigation bar and selecting a patient from the list to modify their programming.

 **Maintenance On**

The patient will start the Exercise by Clicking on "Begin" button.


Begin

Your patient will not wear the Red/Blue glasses for this procedure.


MFBF mode may be selected and your patient will be instructed to wear a Red Eye Patch, which may be ordered from HTS.

They may see

Do not wear your Red/Blue Glasses



Wear the Red Filtering Eye Patch




While running Pursuits

The E will move smoothly across the screen and randomly change its orientation. The task is to indicate when the E changes orientation and is now facing a different direction, either up, down, left, or right.

When the orientation of the E changes, the patient should press the arrow keys on their keyboard in the direction the E is facing. They should respond as rapidly as possible.

Correct responses will be denoted by a BEEP tone. Incorrect responses will be denoted by a BOOP tone.




Percent Correct

A correct response denotes responding to the direction the E is facing before the next orientation change.

The default goal is 80%.

You may move your mouse pointer over the columns in the graph to display the numerical values.

Percent Correct



Goal



Last Session



This Session


Average Response Time

Average Response Time to orientation changes.

The default goal is .95 of a second.

You may move your mouse pointer over the columns in the graph to display the numerical values.

Average Response Time (In Seconds)



Goal



Last Session



This Session

Stars Achieved

After achieving a star, the target speed will increase. After every 4 Stars achieved, the speed will reset and the target size will decrease.

Stars Achieved



Settings

Image Size

Size of Image.

Image Size

1	2	3
---	---	---

MFBF

Monocular In Binocular Field

☒ **MFBF**

**Auditory
Distractors**

Auditory Distractors will present the patient with various verbal prompts that they should verbally respond to. For example, repeat the following word. The patient may select, from a list, the distractor task they wish to respond to.

☒ **Auditory Distractors**

Duration

Assigned Exercise time in minutes.

Duration

-	3	+
---	---	---

Session #

The number of sessions completed.

Session: #1

Viewing Distance

Working distance from the patients eyes to the screen in inches.



