

HTS2

HTS2 - Online VT

Saccades Doctor Manual

© Copyright HTS Inc. 2019



HTS2 - Online VT

Developed By

Dr. Jeffery Cooper, OD, FCOVD
Rodney K. Bortel
Ted L. Conner

Saccades

When in the default Daily Therapy Protocol, Saccades will run for 3 minutes each session.

TODAY'S ASSIGNMENT	MIN.
Pursuits	03:00
Saccades	03:00
Divergence	07:00
Convergence	07:00
Accommodative Rock	05:00

If all of the Stars have been earned for Saccades and you wish your patient to continue Saccades, you may place the Exercise in Maintenance Mode. A Exercise in Maintenance Mode will always be displayed in the patients Daily Therapy Protocol until removed from Maintenance. Maintenance Mode may be accessed by Clicking Existing Patients from the Navigation bar and selecting a patient from the list to modify their programming.

 **Maintenance On**

The patient will start the Exercise by Clicking on "Begin" button.


Begin

Your patient will not wear the Red/Blue glasses for this procedure.


MFBF mode may be selected and your patient will be instructed to wear a Red Eye Patch, which may be ordered from HTS.

They may see

Do not wear your Red/Blue Glasses



Wear the Red Filtering Eye Patch




While running Saccades

The task is to indicate the direction the arrow on the screen is pointing. The patient responds with the left and right arrow keys on their keyboard. They should respond as accurately and rapidly as possible.

After each response, the arrow will randomly jump to a different location on the screen. The goal of this Exercise is to improve both speed and accuracy.

Correct responses will be denoted by a BEEP tone. Incorrect responses will be denoted by a BOOP tone.




Percent Correct

A correct response denotes responding to the direction the arrow is pointing.

The default goal is 80%.

You may move your mouse pointer over the columns in the graph to display the numerical values.

Percent Correct



Goal



Last Session



This Session


Average Response Time

Average Response Time to indicate the direction the arrow is pointing.

The default goal is .95 of a second.

You may move your mouse pointer over the columns in the graph to display the numerical values.

Average Response Time (In Seconds)



Goal



Last Session



This Session

Stars Achieved

After every 4 Stars achieved, the arrow will decrease in size.

Stars Achieved



Settings

Image Size

Size of Arrow.

Image Size

1	2	3
---	---	---

Image Speed

Speed of Image.

Image Speed

1	2	3	4
---	---	---	---

MFBF

Monocular In Binocular Field

☒ MFBF

Auditory Distractors

Auditory Distractors will present the patient with various verbal prompts that they should verbally respond to. For example, repeat the following word. The patient may select, from a list, the distractor task they wish to respond to.

☒ Auditory Distractors

Duration

Assigned Exercise time in minutes.

Duration

-	3	+
---	---	---

Session #

The number of sessions completed.

Session: #1

Viewing Distance

Working distance from the patients eyes to the screen in inches.



