

# HTS2

## HTS2 - Online VT

---

### Vergence Base Down Doctor Manual

© Copyright HTS Inc. 2019



**HTS2 - Online VT**

**Developed By**

Dr. Jeffery Cooper, OD, FCOVD  
Rodney K. Bortel  
Ted L. Conner

## Vergence Base Down

Vergence Base Down may be added to the Daily Therapy Protocol and will run for 7 minutes each session.

TODAY'S ASSIGNMENT	MIN.
Vergence Base Up	07:00
Vergence Base Down	07:00

The Exercise will default to the Maintenance Mode, so no stars will be assigned and it will run until removed from the Daily Protocol.

 **Vergence Base Up Assigned**

The patient will start the Exercise by Clicking on "Begin" button.

**Begin**

## While running Vergence Base Down

Vergence Base Down does not have assigned goals, therefore no stars will be earned.


The Spaceship and Clicker modes are not available for the Vergence Base Down Exercise.

### Classic Mode

The task is, using the arrow keys on the keyboard, to respond to the position of the small square that is popping out of the screen.

The small square will appear either left, right, above, or below the center of the screen.

Correct responses will be denoted by a BEEP tone. Incorrect responses will be denoted by a BOOP tone.




## Results

### Prism Diopter

Vergence Demand measured in Prism Diopters.

You may move your mouse pointer over the columns in the graph to display the numerical values.



# Settings

## Target Size

Size of image, Large or Small.

### Target Size

Large	Small
-------	-------

## Response Time Out

Time allowed to make a response before an error is scored.

### Response Timeout Seconds

-	10	+
---	----	---

## Auditory Distractors

Auditory Distractors will present the patient with various verbal prompts that they should verbally respond to. For example, repeat the following word. The patient may select, from a list, the distractor task they wish to respond to.

### Auditory Distractors

**Duration**

Assigned Exercise time in minutes.

**Duration**

-	7.0	+
---	-----	---

**Session #**

The number of sessions completed.

**Session: #1**

**Viewing Distance**

Working distance from the patients eyes to the screen in inches.



